




# Zephyr<sup>®</sup> Stretchers

Patient Positioning  
and Transfer System



## Zephyr Patient Transport Stretchers

Conventional and MR conditional stretchers are available that are specifically designed for the Zephyr Patient Positioning and Transfer Systems.



## Benefits of the Zephyr Patient Transport Stretcher

- Continuous flat surface with translating top enables smooth patient transfer from imaging/treatment to stretcher
- Structural reinforcements stabilize stretcher during transfers
- Patients may be treated directly on stretcher for HDR Brachytherapy
- Safety stops at head/foot stretcher surface
- Four-wheel braking system

## Benefits of the Zephyr MR Conditional Stretcher

- MR Conditional for 3T
- Continuous flat surface with folding bridge enables smooth patient transfer from imaging/treatment to stretcher
- Designed for stable transfers
- Patients may be treated directly on stretcher for HDR Brachytherapy
- Four-wheel braking system



## Zephyr Stretcher Matrix

Model	ZTG-2	ZTG-3
<b>Description</b>	Zephyr MR Conditional Stretcher	Zephyr Patient Transport Stretcher
<b>Product Focus</b>	Designed for the Zephyr HDR or Zephyr XL Patient Transfer System	Designed for the Zephyr HDR or Zephyr XL Patient Transfer System
<b>Dimensions</b>		
Weight	250 lb (113.6 kg)	340 lb (154.5 kg)
Length	83.25 in (211.5 cm)	83 in (210.8 cm)
Height	32 in (81.3 cm)	MAX: 35.5 in (90.2 cm) MIN: 22.2 in (56.4 cm)
Overall Width (Side rails stored)	24.9 in (63.2 cm)	29.95 in (76.1 cm)
Overall Width (Side rails Up)	35 in (88.9 cm)	32 in (81.3 cm)
Side Rail Length	48 in (121.9 cm)	47 in (119.4 cm)
Wheel Base	24 in x 66 in (61 cm x 167.6 cm)	24 in x 50.5 in (61 cm x 128.3 cm)
Caster Size	5 in (12.7 cm)	8 in (20.3 cm)
Maximum Working Load	600 lb (272 kg)	600 lb (272 kg)
<b>Key Features</b>		
<b>Continuous Flat Surface</b>	Yes	Yes
<b>Bridging Capabilities</b>	Bridge Folds Down Bridge Width: 6 in (15.24 cm)	Translates +/- 4.5 inches (+/- 11.4 cm)
<b>End Stops</b>	Yes	Yes
<b>Height Stops</b>	Fixed Height	Yes
<b>Adjustable Height</b>	Fixed Height	Yes
<b>Trendelenburg / Reverse Trendelenburg</b>	N/A	18°
<b>MR Conditional</b>	Yes	No
<b>SingleAction Locking Mechanism</b>	Yes	Yes
<b>5th wheel provides ease of steering</b>	Yes	Yes
<b>Patient Side Rails</b>	Yes	Yes
<b>IV Pole</b>	No	Yes

**DIACOR** | [www.diacorinc.com](http://www.diacorinc.com)  
 2550 Decker Lake Blvd. Suite 26  
 Salt Lake City, UT 84119  
 Toll free (800) 342.2679  
 Phone (801) 467.0050  
 Fax (801) 487.3258



Precision comes standard.  
 Since 1983.



Check out our latest videos on the web!