




Zephyr[®] Stretchers

Patient Positioning
and Transfer System



Zephyr Patient Transport Stretchers

Conventional and MR conditional stretchers are available that are specifically designed for the Zephyr Patient Positioning and Transfer Systems.



Benefits of the Zephyr Patient Transport Stretcher

- Continuous flat surface with translating top enables smooth patient transfer from imaging/treatment to stretcher
- Structural reinforcements stabilize stretcher during transfers
- Patients may be treated directly on stretcher for HDR Brachytherapy
- Safety stops at head/foot stretcher surface
- Four-wheel braking system

Benefits of the Zephyr MR Conditional Stretcher

- MR Conditional for 3T
- Continuous flat surface with folding bridge enables smooth patient transfer from imaging/treatment to stretcher
- Designed for stable transfers
- Patients may be treated directly on stretcher for HDR Brachytherapy
- Four-wheel braking system



Zephyr Stretcher Matrix

Model	ZTG-2	ZTG-3
Description	Zephyr MR Conditional Stretcher	Zephyr Patient Transport Stretcher
Product Focus	Designed for the Zephyr HDR or Zephyr XL Patient Transfer System	Designed for the Zephyr HDR or Zephyr XL Patient Transfer System
Dimensions		
Weight	250 lb (113.6 kg)	340 lb (154.5 kg)
Length	83.25 in (211.5 cm)	83 in (210.8 cm)
Height	32 in (81.3 cm)	MAX: 35.5 in (90.2 cm) MIN: 22.2 in (56.4 cm)
Overall Width (Side rails stored)	24.9 in (63.2 cm)	29.95 in (76.1 cm)
Overall Width (Side rails Up)	35 in (88.9 cm)	32 in (81.3 cm)
Side Rail Length	48 in (121.9 cm)	47 in (119.4 cm)
Wheel Base	24 in x 66 in (61 cm x 167.6 cm)	24 in x 50.5 in (61 cm x 128.3 cm)
Caster Size	5 in (12.7 cm)	8 in (20.3 cm)
Maximum Working Load	600 lb (272 kg)	600 lb (272 kg)
Key Features		
Continuous Flat Surface	Yes	Yes
Bridging Capabilities	Bridge Folds Down Bridge Width: 6 in (15.24 cm)	Translates +/- 4.5 inches (+/- 11.4 cm)
End Stops	Yes	Yes
Height Stops	Fixed Height	Yes
Adjustable Height	Fixed Height	Yes
Trendelenburg / Reverse Trendelenburg	N/A	18°
MR Conditional	Yes	No
SingleAction Locking Mechanism	Yes	Yes
5th wheel provides ease of steering	Yes	Yes
Patient Side Rails	Yes	Yes
IV Pole	No	Yes

DIACOR | www.diacorinc.com
 2550 Decker Lake Blvd. Suite 26
 Salt Lake City, UT 84119
 Toll free (800) 342.2679
 Phone (801) 467.0050
 Fax (801) 487.3258



Precision comes standard.
 Since 1983.



Check out our latest videos on the web!