




# Zephyr<sup>®</sup> Stretchers

Patient Positioning  
and Transfer System



## Zephyr Patient Transport Stretchers

Conventional and MR conditional stretchers are available that are specifically designed for the Zephyr Patient Positioning and Transfer Systems.



## Benefits of the Zephyr Patient Transport Stretcher

- Continuous flat surface with translating top enables smooth patient transfer from imaging/treatment to stretcher
- Structural reinforcements stabilize stretcher during transfers
- Patients may be treated directly on stretcher for HDR Brachytherapy
- Safety stops at head/foot stretcher surface
- Four-wheel braking system

## Benefits of the Zephyr MR Conditional Stretcher

- MR Conditional for 3T
- Continuous flat surface with folding bridge enables smooth patient transfer from imaging/treatment to stretcher
- Designed for stable transfers
- Patients may be treated directly on stretcher for HDR Brachytherapy
- Four-wheel braking system



## Zephyr Stretcher Matrix

| Model  | ZTG-2  | ZTG-3  |
|--|--|--|
| <b>Description</b>                           | Zephyr MR Conditional Stretcher                                  | Zephyr Patient Transport Stretcher                               |
| <b>Product Focus</b>                         | Designed for the Zephyr HDR or Zephyr XL Patient Transfer System | Designed for the Zephyr HDR or Zephyr XL Patient Transfer System |
| <b>Dimensions</b>                            |  |  |
| Weight                                       | 250 lb (113.6 kg)  | 340 lb (154.5 kg)  |
| Length                                       | 83.25 in (211.5 cm)  | 83 in (210.8 cm)   |
| Height                                       | 32 in (81.3 cm)  | MAX: 35.5 in (90.2 cm)<br>MIN: 22.2 in (56.4 cm)                 |
| Overall Width (Side rails stored)            | 24.9 in (63.2 cm)  | 29.95 in (76.1 cm)   |
| Overall Width (Side rails Up)                | 35 in (88.9 cm)  | 32 in (81.3 cm)  |
| Side Rail Length                             | 48 in (121.9 cm)   | 47 in (119.4 cm)   |
| Wheel Base                                   | 24 in x 66 in<br>(61 cm x 167.6 cm)                              | 24 in x 50.5 in<br>(61 cm x 128.3 cm)                            |
| Caster Size                                  | 5 in (12.7 cm)   | 8 in (20.3 cm)   |
| Maximum Working Load                         | 600 lb (272 kg)  | 600 lb (272 kg)  |
| <b>Key Features</b>                          |  |  |
| <b>Continuous Flat Surface</b>               | Yes  | Yes  |
| <b>Bridging Capabilities</b>                 | Bridge Folds Down<br>Bridge Width: 6 in (15.24 cm)               | Translates +/- 4.5 inches<br>(+/- 11.4 cm)                       |
| <b>End Stops</b>                             | Yes  | Yes  |
| <b>Height Stops</b>                          | Fixed Height   | Yes  |
| <b>Adjustable Height</b>                     | Fixed Height   | Yes  |
| <b>Trendelenburg / Reverse Trendelenburg</b> | N/A  | 18°  |
| <b>MR Conditional</b>                        | Yes  | No   |
| <b>SingleAction Locking Mechanism</b>        | Yes  | Yes  |
| <b>5th wheel provides ease of steering</b>   | Yes  | Yes  |
| <b>Patient Side Rails</b>                    | Yes  | Yes  |
| <b>IV Pole</b>                               | No   | Yes  |

**DIACOR** | [www.diacorinc.com](http://www.diacorinc.com)  
 2550 Decker Lake Blvd. Suite 26  
 Salt Lake City, UT 84119  
 Toll free (800) 342.2679  
 Phone (801) 467.0050  
 Fax (801) 487.3258



Precision comes standard.  
 Since 1983.



Check out our latest videos on the web!